

## **Maximo Application Suite Instrument Standard Transaction in Support of Life Sciences**

### Background

Out of box Maximo Application Suite supports reporting instrument standards (i.e. Tools and corresponding Rotating Assets) on calibration work orders as tools transactions. However, Maximo does not allow to edit such transactions even if those have been reported in error after the approval of the transaction which is either on a save action or work order Complete status change. This negatively impacts instrument standard reverse traceability (i.e. Tool Usage and related electronic reports) and historical work order review.

### User Story 1

Metrology users would benefit from an ability to edit tool transactions (TOOLTRANS object) even after Complete or Closed work order status.

### Solution Requirements

1. Ability to edit work order tool transactions on work orders in Complete status in Maximo Application Suite
2. Ability to edit work order tool transactions on work orders in Closed status (ex. via Edit History menu option) in Maximo Application Suite
3. Ability to control security access to work order application which grants an ability to edit work order tool transactions on work orders in Complete or Closed status in Maximo Application Suite